

foottopia



Foot Massager

USER MANUAL

www.reathlete.com
support@reathlete.com



Thank you for purchasing FOOTTOPIA!

This foot massager combines multiple functions, like pressure, heat, and rubbing, to provide comprehensive relaxation and recovery for your feet. FOOTTOPIA will make you feel better day by day by relieving and alleviating the following issues:

Foot and Toe Pain

Tense Muscles

Cold and Weak Feet

Pain in Tendons and Joints

Foot Fatigue and Numbness

Poor Circulation

Stiffness in The Feet

General Stress and Sleep Problems

Please read this manual before using FOOTTOPIA and keep the manual for reference.

If you have any questions about the foot massager or face any issues when using it, feel free to contact us directly for answers and assistance.

Don't use FOOTTOPIA in any ways that are not described in this manual.

Control Panel

The control panel allows you to conveniently regulate your massage as it happens and to stay in a comfortable position.



Power Button

Forward / Reverse Button

Intensity Button

Heat Button



Plug It In!

Don't forget to power the massager using its cord before you attempt to use the device for your massage session.



How to Use FOOTTOPIA

FOOTTOPIA is built ergonomically to make your feet feel perfectly comfortable. It works different regions of your feet to relieve pain and improve blood flow within minutes.

1. Connect the cord to the massager and plug it into the socket.
2. Put your feet or calves onto the massager.
3. Press the Power button to turn the massager on: it will start in the automatic mode.
4. Press the Heat button to turn the heat on. Press the Heat button again to turn the heat off.
5. Press the Forward / Reverse button to change the spin direction.
6. Press the Intensity button to change the speed of massage. Three intensities are available. Each time you press the Intensity switch, the speed goes up. When the maximum speed is reached, press the Intensity button to go back to the initial speed.
7. Play with the configurations, and build a foot massage routine that will keep your feet healthy and fit every day
8. To turn the device off, press the Power button and unplug the massager.

Troubleshooting

Problem	Cause	Solution
The device doesn't work.	No power. The device hasn't been plugged into the socket. The Power switch hasn't been pressed.	Make sure the device is properly powered and turn it on by pressing the Power switch.
The operation is slow and weak.	The voltage is low.	Use a regulator.
The device is making noise.	The device is new.	After 20 hours of operation, the noise will be decreased.
The device stops working suddenly while massage is in progress.	The device is not properly powered. The Power switch has been accidentally pressed. The overheat protection device was activated.	Make sure the device is plugged in and powered. Press the Power switch to turn the device on. Turn the device off, unplug it, and let it rest and cool down for 20 minutes. If the problem persists, stop using the device: contact our support service.

Safety Precautions

- Do not use any power supply except 110V-220V.
- Unplug the device when you're not using it for massage.
- Stop using the device if it causes discomfort or pain.
- Do not use the device under pillows, blankets, or other covers to avoid fire, electric shock, and injuries.
- If the device is used by or close to children or persons with physical or mental disabilities, they need to be supervised by adult guardians.
- Do not use parts not recommended by the manufacturer on this device. Do not use the device in any ways not described in this manual.
- Do not pull or drag the device by its cord, and do not use the cord as a handle.
- Keep the cord at a safe distance from heated surfaces.
- Do not insert or drop any objects in the openings on the device.
- Do not use the device outdoors.
- Do not use the device in environments where there is a lot of spraying materials or oxygen.
- If the power supply is cut off, turn the device off and unplug it.
- If any part of the device, including the cord and the plug, has been damaged due to falling, contact with water, or other accidents, return the device to the maintenance service center to have the device inspected and repaired.
- Do not use the device in environments with excessive humidity or dense dust to avoid electric shock and failures.
- Do not use the device in environments where the temperature is more than 40° C (104° F).
- Do not place the device near heat sources or in direct sunlight.
- The device is to be used indoors at room temperature. If you move the device from a cold place to a warm

place, its operation can be negatively affected. After you move the device, it is recommended to let it rest at room temperature for one hour before you start using it.

- Do not use the device if its cloth cover is damaged.
- Do not use the device in combination with other curative equipment like electric blankets.
- Avoid using the device within one hour after eating.
- Do not use the device if your hands, feet, or other body parts are wet.
- Continuous operation of the device on the same body part can cause overstimulation and adverse effects. Do not use the device if it causes discomfort.
- This device should not be used for self-diagnosing or treatment. Diagnosis and treatment should be handled by professional medical personnel.
- If you take medications or otherwise receive medical treatment, consult your physician before using the device.
- If you are pregnant, consult your physician before using the device.
- The recommended time limit is 15 minutes. After a massage session, let the device rest for at least 10 minutes. If the device is used continuously for a long period, the overheat protection may be activated, so the device will automatically shut down.
- If the surface of the device is hot, users sensitive to heat need to be supervised when using the device.
- Do not use the device when intoxicated.
- Do not use the device on your head, elbows, knees, chest, back, or abdomen.
- If you spill some liquid on the device, immediately turn it off and unplug it.
- Stop using the device immediately in cases of pain, any mechanical damage, and lightning.

Cleaning and Maintenance

- Turn off and unplug the device before cleaning it.
- Clean the device with a soft cloth.
- Do not use erosive detergents. You can use a neutral detergent to clean the device.
- Do not scratch any part of the device with sharp objects.
- Do not put the device in a wet environment.
- Keep the device in places with good ventilation.
- Do not exert excessive pressure on the device.

Technical Details

Rated voltage: 12 V

Rated power: 24 W

Suggested massage time: 15 min

Executive standard:

GB4706.1-2005, GB4706.10-2008

foottopia

**If you encounter
any other problems,
please contact us at
support@reathlete.com**